



perform[®] sterile cleaner ND

Sterile, ready to use solution based on non-ionic surfactants for the cleaning of surfaces in cleanrooms



- Sterile (γ -irradiated) and microbial filtrated (0,2 μ m)
- Suitable for cleanrooms grade A/B, ISO 5/6
- Ready to use in a handy, sterile 500ml spray system
- Color and perfume free
- Excellent cleaning performance

Application areas

Disinfection of surfaces and equipment / devices in cleanrooms, especially for class A and B cleanroom areas, according to:

- Annex 1 of EU Guide to Good Manufacturing practice: Manufacture of sterile medicinal products, March 2009.
- FDA Guidance for Industry: Sterile drug products produced by aseptic processing – current Good Manufacturing practice; U.S. Department of Health and Human Services, Food and Drug Administration, Sept. 2004

Product data

Labelling according to Regulation (EC) No. 648/2004: <5% non-ionic surfactants

Chemical-physical data:

Density	: appr. 1,00 g/cm ³ / 20°C
Color	: Colorless
pH value	: appr. 5,8 / 20°C / concentrate
Flash point	: Not applicable
Form	: Liquid

Special advice

The sterile 500ml spray system offers the following advantages:

- Convenient size for handling and storage
- Ergonomically shaped spray pump for ease of use
- Remains sterile after opening

The bottles are sterile (with irradiation indicator) and double wrapped individually.

Instructions for use

Apply perform[®] sterile cleaner ND undiluted to objects and surfaces and wipe if necessary.

Packaging unit : 10 x 500 ml bottle
Art.-no. : RREIN-PER-STND

hans j. michael gmbh

Gewerbegebiet Hart 11, D-71554 Weissach i.T., Phone: 0049-7191-9105-0, Fax: 0049-7191-9105-19
e-mail: office@hjm-reinraum.de, website: <http://www.cleanroom-hjm.com>