



## perform<sup>®</sup> classic concentrate Q-Plus

Liquid concentrate based on a combination of quaternary ammonium compounds, alcohols, amphoteric glycine derivatives and surfactants for disinfecting and cleaning surfaces



- Usable in cleanroom C/D (ISO 7/8)
  - Fast and effective
- Excellent material compatibility
  - Colour and perfume free
- Efficacy tested according to European norms (EN)

### **Application areas**

For disinfection and cleaning of all types of surfaces. Suitable for use in all areas where, in addition to a broad spectrum of efficacy, strict requirements are placed on the environmental behaviour and user friendliness of a surface disinfectant.

### **Composition**

100g solution contain the following active ingredients:

22g Alkyl(C12-16)dimethylbenzylammoniumchlorid (ADBAC/BKC (C12-16)), 17g 2-Phenoxyethanol, 0,9g Amines, N-C12-14-(even-numbered)-alkyltrimethylenedi-, reaction products with chloroacetic acid.

Labelling according to Regulation (EC) No. 648/2004: 5-15% non-ionic surfactants

### **Microbiological efficacy**

<b>Efficacy:</b>	<b>Concentration:</b>	<b>Contact time:</b>
<b>Bactericidal, EN1276, EN13697</b>	0,5% (5ml/l)	60 minutes
<b>Fungicidal, EN1650, EN13697</b>		
<b>Enveloped viruses (incl. HIV, HBV and HCV), in acc. with DVV</b>		

### **Chemical-physical data**

Density	: Appr. 1,0 g/cm <sup>3</sup> / 20°C
Colour	: Nearly colourless
pH-value	: Appr. 8,6 / 20°C (Concentrate)
Flashpoint	: 51°C / Method: DIN 51755 Part 1
Form	: Liquid
Viscosity, dynamic	: Appr. 21 mPa*s / 20°C / Method: ISO 3219

### **Special advice**

Use biocides safely. Always read the label and product information before use. Wearing gloves (e.g. made of butyl rubber) is recommended.

IE: PCS number: 95502



### **Instructions for use**

Heavily contaminated surfaces are to be pre-cleaned. Prepare the working solution by filling water of an appropriate quality into a tub or bucket. Add the concentrate to the water and stir gently. Wipe surfaces, floors and other wipeable areas with a moist cloth.

perform® classic concentrate Q-Plus has excellent cleaning properties, therefore the use of additional cleaning additives is unnecessary. Depending on the type and quality of the floor, supplementary floor care may be necessary. Do not combine this product with other products. This may lead to reactions with anionic surfactants due to the content of cationic active ingredients.

<b>Art.-no.:</b>	<b>PU:</b>
RREIN-PER-CLQP	5 Liters canister

**hans j. michael gmbh**

Gewerbegebiet Hart 11, D-71554 Weissach i.T., Phone: 0049-7191-9105-0, Fax: 0049-7191-9105-19  
e-mail: [office@hjm-reinraum.de](mailto:office@hjm-reinraum.de), website: <http://www.cleanroom-hjm.com>